

Items required for each test		Corrosion		Oxidation	
		ASTM D130	ASTM D4048	ASTM D525, D873	ASTM D5304
		IP 154	IP 112	IP 40, IP 138	
11241-0	Silicone Carbide Paper 150 grit (pack of 50)	●	●		
11420-0	Support Assembly (pair)		●		
11460-0	Silicon Carbide Paper, 65 micron, P240 FEPA Grade (pack of 50)	●	●		
11470-0	Silicon Carbide Paper P100 FEPA Grade (pack of 50)	●	●		
11480-0	Silicon Carbide Grains 150 Mesh 500 g (pack of 5)	●	●		
11500-0	Pressure Vessel	●			
11514-0	Tongs	●	●		
11530-0	Multi-Strip Vice	●	●		
11550-0	Copper Test Strip (pack of 30)	●	●		
11570-0	Flat Glass Test Tube (pack of 10)	●	●		
11580-0	ASTM Copper Strip Corrosion Standard	●	●		
11590-0	Test Tube (pack of 10)		●		
ASTMS12C	Thermometer	●	●		
15600-0	Seta Pressure Gauge			▲	▲
15802-2	Oxidation Vessel			●	
15803-0	Vessel Bench Clamp and Holder			●	●
15805-0	Spanner			●	●
15810-0	Liner with Cover (pack of 3)			●	●
99100-2	Analytical Balance			●	●
99200-3	Laboratory Drying Oven			●	●
ASTMS22C	Low Haz Thermometer with Certificate of Conformity			●	
ASTMST61C	Low Haz Thermometer with Works Certificate				●

Items marked ▲ are options for pressure monitoring